

产品名称: AGR3 Mouse Monoclonal Antibody  
产品货号: AMM83147



## 产品概述 (Summary)

|                        |                                |
|------------------------|--------------------------------|
| 产品名称 (Production Name) | AGR3 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody      |
| 宿主 (Host)              | Mouse                          |
| 应用 (Application)       | WB,IHC                         |
| 种属反应性 (Reactivity)     | Human                          |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                    |
|------------------------|----------------------------------------------------|
| 基因名 (Gene Name)        | AGR3                                               |
| 别名 (Alternative Names) | AGR3                                               |
| 基因 ID (Gene ID)        |                                                    |
| 蛋白 ID (SwissProt ID)   | Q8TD06.Purified recombinant fragment of human AGR3 |

## 产品应用 (Application)

|                          |                                  |
|--------------------------|----------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:500-1:1000 |
| 蛋白分子量 (Molecular Weight) | 19.2kDa                          |


## 研究背景 (Background)

产品名称: AGR3 Mouse Monoclonal Antibody  
产品货号: AMM83147




## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis using AGR3 mouse mAb against T-47D lysate



IHC-P analysis of fallopian tube tissue by AGR3 antibody (AMM83147)

## 注意事项 (Note)

For research use only .