

产品名称: TGM2 Mouse Monoclonal Antibody
产品货号: AMM83158



产品概述 (Summary)

产品名称 (Production Name)	TGM2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TGM2
别名 (Alternative Names)	TGM2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P21980.Purified recombinant fragment of human TGM2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	90kDa


研究背景 (Background)

产品名称: TGM2 Mouse Monoclonal Antibody
产品货号: AMM83158




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using TGM2 mouse mAb against A549 lysate



Lane 1: A549 lysate ; Lane 2: TGM2 immunoprecipitated from A549 lysate by AMM83158 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83158 can immunoprecipitate TGM2;



IHC-P analysis of colon tissue by TGM2 antibody (AMM83158)

注意事项 (Note)

产品名称: TGM2 Mouse Monoclonal Antibody
产品货号: AMM83158



For research use only .