

产品名称: ANXA3 Mouse Monoclonal Antibody
产品货号: AMM83159



产品概述 (Summary)

产品名称 (Production Name)	ANXA3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA3
别名 (Alternative Names)	ANXA3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P12429.Purified recombinant fragment of human ANXA3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	36kDa


研究背景 (Background)

产品名称: ANXA3 Mouse Monoclonal Antibody
产品货号: AMM83159




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using ANXA3 mouse mAb against HeLa lysate



Lane 1: HeLa lysate ; Lane 2: ANXA3 immunoprecipitated from HeLa lysate by AMM83159 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83159 can immunoprecipitate ANXA3;



IHC-P analysis of kidney tissue by ANXA3 antibody (AMM83159)

注意事项 (Note)

产品名称: ANXA3 Mouse Monoclonal Antibody
产品货号: AMM83159



For research use only .