

产品名称: TRIM34 Mouse Monoclonal Antibody
产品货号: AMM83167



产品概述 (Summary)

产品名称 (Production Name)	TRIM34 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIM34
别名 (Alternative Names)	TRIM34
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BYJ4.Purified recombinant fragment of human TRIM34

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	56.9kDa


研究背景 (Background)

产品名称: TRIM34 Mouse Monoclonal Antibody
产品货号: AMM83167




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using TRIM34 mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: TRIM34 immunoprecipitated from Jurkat lysate by AMM83167 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83167 can immunoprecipitate TRIM34;



IHC-P analysis of placenta tissue by TRIM34 antibody (AMM83167)

注意事项 (Note)

产品名称: TRIM34 Mouse Monoclonal Antibody
产品货号: AMM83167



For research use only .