

产品名称: VEGFC Mouse Monoclonal Antibody
产品货号: AMM83177



产品概述 (Summary)

产品名称 (Production Name)	VEGFC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VEGFC
别名 (Alternative Names)	VEGFC
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49767.Purified recombinant fragment of human VEGFC

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 42 kDa

研究背景 (Background)

产品名称: VEGFC Mouse Monoclonal Antibody
产品货号: AMM83177




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using VEGFC mouse mAb against rat heart tissue lysate



Lane 1: Rat heart lysate ; Lane 2: VEGFC immunoprecipitated from rat heart lysate by AMM83177 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83177 can immunoprecipitate VEGFC;

注意事项 (Note)

For research use only .