

产品名称: PTGES2 Mouse Monoclonal Antibody
产品货号: AMM83185



产品概述 (Summary)

产品名称 (Production Name)	PTGES2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTGES2
别名 (Alternative Names)	PTGES2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H7Z7.Purified recombinant fragment of human PTGES2

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: PTGES2 Mouse Monoclonal Antibody
产品货号: AMM83185




研究领域 (Research Area)

图片 (Image Data)



Lane 1: HeLa lysate ; Lane 2: PTGES2 immunoprecipitated from HeLa lysate by AMM83185 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83185 can immunoprecipitate PTGES2;



IHC-P analysis of heart muscle tissue by PTGES2 antibody (AMM83185)

注意事项 (Note)

For research use only .