

产品名称: UBE2Z Mouse Monoclonal Antibody
产品货号: AMM83213



产品概述 (Summary)

产品名称 (Production Name)	UBE2Z Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE2Z
别名 (Alternative Names)	UBE2Z
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H832.Purified recombinant fragment of human UBE2Z

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	38.2kDa


研究背景 (Background)

产品名称: UBE2Z Mouse Monoclonal Antibody
产品货号: AMM83213




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using UBE2Z mouse mAb against MCF-7 lysate



Lane 1: MCF7 lysate ; Lane 2: UBE2Z immunoprecipitated from MCF7 lysate by AMM83213 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83213 can immunoprecipitate UBE2Z;



IHC-P analysis of testis tissue by UBE2Z antibody (AMM83213)

注意事项 (Note)

产品名称: UBE2Z Mouse Monoclonal Antibody
产品货号: AMM83213



For research use only .