

产品名称: MTHFD2 Mouse Monoclonal Antibody
产品货号: AMM83224



产品概述 (Summary)

产品名称 (Production Name)	MTHFD2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MTHFD2
别名 (Alternative Names)	MTHFD2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P13995.Purified recombinant fragment of human MTHFD2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	37.9kDa


研究背景 (Background)

产品名称: MTHFD2 Mouse Monoclonal Antibody
产品货号: AMM83224




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using MTHFD2 mouse mAb against A-431 lysate



Lane 1: A-431 lysate ; Lane 2: MTHFD2 immunoprecipitated from A-431 lysate by AMM83224 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83224 can immunoprecipitate MTHFD2;



IHC-P analysis of kidney tissue by MTHFD2 antibody (AMM83224)

注意事项 (Note)

产品名称: MTHFD2 Mouse Monoclonal Antibody
产品货号: AMM83224



For research use only .