

产品名称: SCGB3A2 Mouse Monoclonal Antibody
产品货号: AMM83229



产品概述 (Summary)

产品名称 (Production Name)	SCGB3A2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCGB3A2
别名 (Alternative Names)	SCGB3A2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96PL1.Purified recombinant fragment of human SCGB3A2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	10.2kDa


研究背景 (Background)

产品名称: SCGB3A2 Mouse Monoclonal Antibody
产品货号: AMM83229



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using SCGB3A2 mouse mAb against mouse lung tissue lysate

注意事项 (Note)

For research use only .