

产品名称: Mitf (phospho S157) Mouse Monoclonal Antibody
产品货号: AMM83238



产品概述 (Summary)

产品名称 (Production Name)	Mitf (phospho S157) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mitf (phospho S157)
别名 (Alternative Names)	Mitf (phospho S157)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O88368.Mitf peptide

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: Mitf (phospho S157) Mouse Monoclonal Antibody
产品货号: AMM83238



研究领域 (Research Area)

图片 (Image Data)



IHC-P analysis of cutaneous melanoma tissue by Mitf antibody (AMM83238)

注意事项 (Note)

For research use only .