

产品名称: CAPS Mouse Monoclonal Antibody
产品货号: AMM83245



产品概述 (Summary)

产品名称 (Production Name)	CAPS Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAPS
别名 (Alternative Names)	CAPS
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13938.Purified recombinant fragment of human CAPS

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	30.2kDa


研究背景 (Background)

产品名称: CAPS Mouse Monoclonal Antibody
产品货号: AMM83245




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using CAPS mouse mAb against U-251MG lysate



Lane 1: U-251MG lysate ; Lane 2: CAPS immunoprecipitated from U-251MG lysate by AMM83245 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83245 can immunoprecipitate CAPS;



IHC-P analysis of fallopian tube tissue by CAPS antibody (AMM83245)

注意事项 (Note)

产品名称: CAPS Mouse Monoclonal Antibody
产品货号: AMM83245



For research use only .