

产品名称: AKR1B1 Mouse Monoclonal Antibody
产品货号: AMM83248



产品概述 (Summary)

产品名称 (Production Name)	AKR1B1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AKR1B1
别名 (Alternative Names)	AKR1B1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P15121.Purified recombinant fragment of human AKR1B1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	36kDa


研究背景 (Background)

产品名称: AKR1B1 Mouse Monoclonal Antibody
产品货号: AMM83248




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using AKR1B1 mouse mAb against HeLa lysate



Lane 1: HeLa lysate ; Lane 2: AKR1B1 immunoprecipitated from HeLa lysate by AMM83248 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83248 can immunoprecipitate AKR1B1;



IHC-P analysis of adrenal gland tissue by AKR1B1 antibody (AMM83248)

注意事项 (Note)

产品名称: AKR1B1 Mouse Monoclonal Antibody
产品货号: AMM83248



For research use only .