

产品名称: Acetyl-SMC2 (Lys114) Mouse Monoclonal Antibody
产品货号: AMM83264



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-SMC2 (Lys114) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMC2 (acetyl K114)
别名 (Alternative Names)	SMC2 (acetyl K114)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95347.Human SMC2 peptide

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:100-1:200
蛋白分子量 (Molecular Weight)	135.7kDa


研究背景 (Background)

产品名称: Acetyl-SMC2 (Lys114) Mouse Monoclonal Antibody
产品货号: AMM83264




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using SMC2 mouse mAb against Raji lysate



Lane 1: Raji lysate; Lane 2: SMC2 immunoprecipitated from Raji lysate by P151963; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: P151963 can immunoprecipitate SMC2;

注意事项 (Note)

For research use only .