

产品名称: BCL2L15 Mouse Monoclonal Antibody
产品货号: AMM83266



产品概述 (Summary)

产品名称 (Production Name)	BCL2L15 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL2L15
别名 (Alternative Names)	BCL2L15
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q5TBC7.Purified recombinant fragment of human BCL2L15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	17.7kDa


研究背景 (Background)

产品名称: BCL2L15 Mouse Monoclonal Antibody
产品货号: AMM83266




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using BCL2L15 mouse mAb against HT-29 lysate



Lane 1: HT-29 lysate ; Lane 2: BCL2L15 immunoprecipitated from HT-29 lysate by AMM83266 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83266 can immunoprecipitate BCL2L15;



IHC-P analysis of colon tissue by BCL2L15 antibody (AMM83266)

注意事项 (Note)

产品名称: BCL2L15 Mouse Monoclonal Antibody
产品货号: AMM83266



For research use only .