

产品名称: UBR7 Mouse Monoclonal Antibody
产品货号: AMM83275



产品概述 (Summary)

产品名称 (Production Name)	UBR7 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBR7
别名 (Alternative Names)	UBR7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8N806.Purified recombinant fragment of human UBR7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	48kDa


研究背景 (Background)

产品名称: UBR7 Mouse Monoclonal Antibody
产品货号: AMM83275




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using UBR7 mouse mAb against T-47D lysate



Lane 1: T-47D lysate ; Lane 2: UBR7 immunoprecipitated from T-47D lysate by AMM83275 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83275 can immunoprecipitate UBR7;



IHC-P analysis of testis tissue by UBR7 antibody (AMM83275)

注意事项 (Note)

产品名称: UBR7 Mouse Monoclonal Antibody
产品货号: AMM83275



For research use only .