

产品名称: DTX3L Mouse Monoclonal Antibody
产品货号: AMM83276



产品概述 (Summary)

产品名称 (Production Name)	DTX3L Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DTX3L
别名 (Alternative Names)	DTX3L
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8TDB6.Purified recombinant fragment of human DTX3L

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	83.6kDa


研究背景 (Background)

产品名称: DTX3L Mouse Monoclonal Antibody
产品货号: AMM83276




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using DTX3L mouse mAb against A-431 lysate



Lane 1: A-431 lysate ; Lane 2: DTX3L immunoprecipitated from A-431 lysate by AMM83276 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83276 can immunoprecipitate DTX3L;

注意事项 (Note)

For research use only .