

产品名称: ACE2 Mouse Monoclonal Antibody
产品货号: AMM83290



产品概述 (Summary)

产品名称 (Production Name)	ACE2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACE2
别名 (Alternative Names)	ACE2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BYF1.Purified recombinant fragment of human ACE2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	92kDa


研究背景 (Background)

产品名称: ACE2 Mouse Monoclonal Antibody
产品货号: AMM83290




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using ACE2 mouse mAb against Calu-3 lysate



IHC-P analysis of kidney tissue by ACE2 antibody (AMM83290)

注意事项 (Note)

For research use only .