

产品名称: CBLL2 Mouse Monoclonal Antibody
产品货号: AMM83291



产品概述 (Summary)

产品名称 (Production Name)	CBLL2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CBLL2
别名 (Alternative Names)	CBLL2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8N7E2.Purified recombinant fragment of human CBLL2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	48.8kDa


研究背景 (Background)

产品名称: CBLL2 Mouse Monoclonal Antibody
产品货号: AMM83291




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using CBLL2 mouse mAb against mouse testis tissue lysate



Lane 1: Mouse testis tissue lysate ; Lane 2: CBLL2 immunoprecipitated from mouse testis tissue lysate by AMM83291 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83291 can immunoprecipitate CBLL2;

注意事项 (Note)

For research use only .