

产品名称: NDEL1 Mouse Monoclonal Antibody
产品货号: AMM83294



产品概述 (Summary)

产品名称 (Production Name)	NDEL1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDEL1
别名 (Alternative Names)	NDEL1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9GZM8.Purified recombinant fragment of human NDEL1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 38 kDa


研究背景 (Background)

产品名称: NDEL1 Mouse Monoclonal Antibody
产品货号: AMM83294




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using NDEL1 mouse mAb against mouse brain and rat brain



Lane 1: Mouse brain tissue lysate ; Lane 2: NDEL1 immunoprecipitated from mouse brain tissue lysate by AMM83294 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83294 can immunoprecipitate NDEL1;

注意事项 (Note)

For research use only .