

产品名称: FGF20 Mouse Monoclonal Antibody
产品货号: AMM83296



产品概述 (Summary)

产品名称 (Production Name)	FGF20 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGF20
别名 (Alternative Names)	FGF20
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NP95.Purified recombinant fragment of human FGF20

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa


研究背景 (Background)

产品名称: FGF20 Mouse Monoclonal Antibody
产品货号: AMM83296



研究领域 (Research Area)

图片 (Image Data)



Result: AMM83296 can detect FGF20 by Western blotting.

注意事项 (Note)

For research use only .