

产品名称: ZNF598 Mouse Monoclonal Antibody
产品货号: AMM83297



产品概述 (Summary)

产品名称 (Production Name)	ZNF598 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF598
别名 (Alternative Names)	ZNF598
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q86UK7.Purified recombinant fragment of human ZNF598

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	99kDa


研究背景 (Background)

产品名称: ZNF598 Mouse Monoclonal Antibody
产品货号: AMM83297




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using ZNF598 mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: ZNF598 immunoprecipitated from Jurkat lysate by AMM83297 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83297 can immunoprecipitate ZNF598;

注意事项 (Note)

For research use only .