

产品名称: DBN1 Mouse Monoclonal Antibody
产品货号: AMM83299



产品概述 (Summary)

产品名称 (Production Name)	DBN1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DBN1
别名 (Alternative Names)	DBN1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q16643.Purified recombinant fragment of human DBN1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	71.4kDa


研究背景 (Background)

产品名称: DBN1 Mouse Monoclonal Antibody
产品货号: AMM83299




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 μ g of SH-SY5Y lysate; Lane 2: 15 μ g of rat brain; Lane 3: 15 μ g of mouse brain; Result: AMM83299 can detect DBN1 by Western blotting;



Lane 1: Mouse brain tissue lysate ; Lane 2: DBN1 immunoprecipitated from mouse brain tissue lysate by AMM83299 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ; Result: AMM83299 can immunoprecipitate DBN1;



IHC-P analysis of brain tissue by DBN1 antibody (AMM83299)

产品名称: DBN1 Mouse Monoclonal Antibody
产品货号: AMM83299



注意事项 (Note)

For research use only .