

产品名称: TOMM20 Mouse Monoclonal Antibody
产品货号: AMM83301



产品概述 (Summary)

产品名称 (Production Name)	TOMM20 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TOMM20
别名 (Alternative Names)	TOMM20
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15388.Purified recombinant fragment of human TOMM20

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa


研究背景 (Background)

产品名称: TOMM20 Mouse Monoclonal Antibody
产品货号: AMM83301




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using TOMM20 mouse mAb against mouse brain and rat brain



Western blot analysis using TOMM20 mouse mAb against HEK-293 lysate



IHC-P analysis of kidney tissue by TOMM20 antibody (AMM83301)

注意事项 (Note)

For research use only .