

**产品名称: SARS-CoV-2 Nucleocapsid Protein Mouse
Monoclonal Antibody**
产品货号: AMM83308



产品概述 (Summary)

产品名称 (Production Name)	SARS-CoV-2 Nucleocapsid Protein Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	SARS-CoV-2

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SARS-CoV-2 Nucleocapsid Protein
别名 (Alternative Names)	SARS-CoV-2 Nucleocapsid Protein
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P0DTC9.Recombinant SARS-CoV-2 nucleocapsid protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	55kDa


研究背景 (Background)

产品名称: SARS-CoV-2 Nucleocapsid Protein Mouse
Monoclonal Antibody
产品货号: AMM83308



研究领域 (Research Area)

图片 (Image Data)



Result: AMM83308 can detect recombinant nucleocapsid protein by Western blotting.

注意事项 (Note)

For research use only .