

产品名称: DPCD Mouse Monoclonal Antibody
产品货号: AMM83315



产品概述 (Summary)

产品名称 (Production Name)	DPCD Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPCD
别名 (Alternative Names)	DPCD
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BVM2.Purified recombinant fragment of human DPCD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	23.2kDa


研究背景 (Background)

产品名称: DPCD Mouse Monoclonal Antibody
产品货号: AMM83315




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using DPCD mouse mAb against rat testis lysate



Lane 1: Mouse testis tissue lysate ; Lane 2: DPCD immunoprecipitated from mouse testis tissue lysate by AMM83315 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83315 can immunoprecipitate DPCD;

注意事项 (Note)

For research use only .