

产品名称: NAGK Mouse Monoclonal Antibody
产品货号: AMM83325



产品概述 (Summary)

产品名称 (Production Name)	NAGK Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NAGK
别名 (Alternative Names)	NAGK
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UJ70.Purified recombinant fragment of human NAGK

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	37.4kDa


研究背景 (Background)

产品名称: NAGK Mouse Monoclonal Antibody
产品货号: AMM83325




研究领域 (Research Area)

图片 (Image Data)




Lane 1: 15 μ g of HeLa lysate; Lane 2: 15 μ g of Jurkat lysate; Lane 3: 15 μ g of HEK-293 lysate; Lane 4: 15 μ g of A-431 lysate;

Result: AMM83325 can detect NAGK by Western blotting.



Lane 1: HeLa lysate; Lane 2: NAGK immunoprecipitated from HeLa lysate by AMM83325; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83325 can immunoprecipitate NAGK;



IHC-P analysis of tonsil tissue by NAGK antibody (AMM83325)

产品名称: NAGK Mouse Monoclonal Antibody
产品货号: AMM83325



注意事项 (Note)

For research use only .