

产品名称: TACO1 Mouse Monoclonal Antibody
产品货号: AMM83331



产品概述 (Summary)

产品名称 (Production Name)	TACO1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TACO1
别名 (Alternative Names)	TACO1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BSH4.Purified recombinant fragment of human TACO1

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: TACO1 Mouse Monoclonal Antibody
产品货号: AMM83331




研究领域 (Research Area)

图片 (Image Data)



Lane 1: U-251MG lysate; Lane 2: TACO1 immunoprecipitated from U-251MG by AMM83331; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83331 can immunoprecipitate TACO1;



IHC-P analysis of kidney tissue by TACO1 antibody (AMM83331)

注意事项 (Note)

For research use only .