

产品名称: BASP1 Mouse Monoclonal Antibody
产品货号: AMM83334



产品概述 (Summary)

产品名称 (Production Name)	BASP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BASP1
别名 (Alternative Names)	BASP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P80723.Purified recombinant fragment of human BASP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	22.7kDa


研究背景 (Background)

产品名称: BASP1 Mouse Monoclonal Antibody
产品货号: AMM83334




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using BASP1 mouse mAb against HeLa lysate



IHC-P analysis of brain tissue by BASP1 antibody (AMM83334)

注意事项 (Note)

For research use only .