

产品名称: NPTX1 Mouse Monoclonal Antibody
产品货号: AMM83338



产品概述 (Summary)

产品名称 (Production Name)	NPTX1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPTX1
别名 (Alternative Names)	NPTX1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15818.Purified recombinant fragment of human NPTX1

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: NPTX1 Mouse Monoclonal Antibody
产品货号: AMM83338




研究领域 (Research Area)

图片 (Image Data)



Lane 1: Rat brain tissue lysate ; Lane 2: NPTX1 immunoprecipitated from rat brain tissue lysate by AMM83338 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83338 can immunoprecipitate NPTX1;



IHC-P analysis of brain tissue by NPTX1 antibody (AMM83338)

注意事项 (Note)

For research use only .