

产品名称: LGALS3 Mouse Monoclonal Antibody
产品货号: AMM83345



产品概述 (Summary)

产品名称 (Production Name)	LGALS3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGALS3
别名 (Alternative Names)	LGALS3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P17931.Purified recombinant fragment of human LGALS3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	33kDa


研究背景 (Background)

产品名称: LGALS3 Mouse Monoclonal Antibody
产品货号: AMM83345




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using LGALS3 mouse mAb against HeLa lysate



Lane 1: HeLa lysate ; Lane 2: LGALS3 immunoprecipitated from HeLa lysate by AMM83345 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83345 can immunoprecipitate LGALS3;

注意事项 (Note)

For research use only .