

产品名称: IFIH1 Mouse Monoclonal Antibody
产品货号: AMM83359



产品概述 (Summary)

产品名称 (Production Name)	IFIH1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFIH1
别名 (Alternative Names)	IFIH1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BYX4.Purified recombinant fragment of human IFIH1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	120kDa


研究背景 (Background)

产品名称: IFIH1 Mouse Monoclonal Antibody
产品货号: AMM83359




研究领域 (Research Area)


图片 (Image Data)



Cell lysate of THP-1 (untreated or treated with LPS) was run on 6-18% SDS-PAGE under reducing conditions and blotted onto nitrocellulose membrane. Result: AMM83359 can detect IFIH1 by Western blotting.



Lane 1: LPS treated THP-1 lysate; Lane 2: IFIH1 immunoprecipitated from LPS treated THP-1 lysate by AMM83359; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83359 can immunoprecipitate IFIH1;



IHC-P analysis of thyroid gland tissue by IFIH1 antibody (AMM83359)

产品名称: IFIH1 Mouse Monoclonal Antibody
产品货号: AMM83359



注意事项 (Note)

For research use only .