

产品名称: FGF12 Mouse Monoclonal Antibody
产品货号: AMM83372



产品概述 (Summary)

产品名称 (Production Name)	FGF12 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGF12
别名 (Alternative Names)	FGF12
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61328.Purified recombinant fragment of human FGF12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	27kDa


研究背景 (Background)

产品名称: FGF12 Mouse Monoclonal Antibody
产品货号: AMM83372




研究领域 (Research Area)

图片 (Image Data)



Lane 1: 15 μ g of rat brain tissue lysate; Lane 2: 15 μ g of mouse brain tissue lysate; Lane 3: 15 μ g of brain tissue lysate; Result: AMM83372 can detect FGF12 by Western blotting;



Lane 1: Rat brain lysate ; Lane 2: FGF12 immunoprecipitated from rat brain lysate by AMM83372 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ; Result: AMM83372 can immunoprecipitate FGF12;

注意事项 (Note)

For research use only .