

产品名称: UPB1 Mouse Monoclonal Antibody
产品货号: AMM83393



产品概述 (Summary)

产品名称 (Production Name)	UPB1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UPB1
别名 (Alternative Names)	UPB1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBR1.Purified recombinant fragment of human UPB1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:;Observed MW:42kD


研究背景 (Background)

产品名称: UPB1 Mouse Monoclonal Antibody
产品货号: AMM83393




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using UPB1 mouse mAb against rat liver tissue lysate



IHC-P analysis of liver tissue by UPB1 antibody (AMM83393)

注意事项 (Note)

For research use only .