

产品名称: ARHGEF1 Mouse Monoclonal Antibody
产品货号: AMM83395



产品概述 (Summary)

产品名称 (Production Name)	ARHGEF1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARHGEF1
别名 (Alternative Names)	ARHGEF1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92888.Purified recombinant fragment of human ARHGEF1

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: ARHGEF1 Mouse Monoclonal Antibody
产品货号: AMM83395




研究领域 (Research Area)

图片 (Image Data)



Lane 1: HEK-293 lysate; Lane 2: ARHGEF1 immunoprecipitated from HEK-293 lysate by AMM83395; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83395 can immunoprecipitate ARHGEF1;



IHC-P analysis of spleen tissue by ARHGEF1 antibody (AMM83395)

注意事项 (Note)

For research use only .