

产品名称: Tau Mouse Monoclonal Antibody
产品货号: AMM83407



产品概述 (Summary)

产品名称 (Production Name)	Tau Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tau
别名 (Alternative Names)	Tau
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P10636-8.Purified recombinant fragment of human Tau

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:250
蛋白分子量 (Molecular Weight)	65kDa


研究背景 (Background)

产品名称: Tau Mouse Monoclonal Antibody
产品货号: AMM83407




研究领域 (Research Area)

图片 (Image Data)




Cross reactivity Microtiter wells were coated with recombinant human Tau and phosphorylated Tau proteins. AMM83407 was used as primary antibody and peroxidase conjugated goat anti- mouse IgG was used as the secondary antibody. Result:

AMM83407 can react with human Tau and phosphorylated Tau proteins.



Result: AMM83407 can detect recombinant human Tau by Western blotting.



IHC-P analysis of cerebral cortex tissue by Tau antibody (AMM83407)

产品名称: Tau Mouse Monoclonal Antibody
产品货号: AMM83407



注意事项 (Note)

For research use only .