

产品名称: TRAIP Mouse Monoclonal Antibody
产品货号: AMM83411



产品概述 (Summary)

产品名称 (Production Name)	TRAIP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAIP
别名 (Alternative Names)	TRAIP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BWF2.Purified recombinant fragment of human TRAIP

产品应用 (Application)

稀释比 (Dilution Ratio)	IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: TRAIP Mouse Monoclonal Antibody
产品货号: AMM83411



研究领域 (Research Area)

图片 (Image Data)



Lane 1: Mouse brain tissue lysate; Lane 2: TRAIP immunoprecipitated from mouse brain tissue lysate by AMM83411; Lane 3:
The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83411 can immunoprecipitate TRAIP;

注意事项 (Note)

For research use only .