

产品名称: VASH1 Mouse Monoclonal Antibody
产品货号: AMM83443



产品概述 (Summary)

产品名称 (Production Name)	VASH1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VASH1
别名 (Alternative Names)	VASH1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q7L8A9.Purified recombinant fragment of human VASH1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	40kDa


研究背景 (Background)

产品名称: VASH1 Mouse Monoclonal Antibody
产品货号: AMM83443



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using VASH1 mouse mAb against Jurkat lysate and 15 μ g

注意事项 (Note)

For research use only .