

产品名称: A2M Mouse Monoclonal Antibody
产品货号: AMM83456



产品概述 (Summary)

产品名称 (Production Name)	A2M Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	A2M
别名 (Alternative Names)	A2M
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01023.Alpha 2 macroglobulin purified from human plasma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	170kDa


研究背景 (Background)

产品名称: A2M Mouse Monoclonal Antibody
产品货号: AMM83456




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using human A2M mouse mAb against Hep G2 lysate human plasma



Lane 1: 15 μ g of Hep G2 lysate; Lane 2: A2M immunoprecipitated from Hep G2 lysate by AMM83456; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83456 can immunoprecipitate A2M;

注意事项 (Note)

For research use only .