

产品名称: DCBLD2 Mouse Monoclonal Antibody
产品货号: AMM83487



产品概述 (Summary)

产品名称 (Production Name)	DCBLD2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DCBLD2
别名 (Alternative Names)	DCBLD2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96PD2.Purified recombinant fragment of human DCBLD2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	85kDa


研究背景 (Background)

产品名称: DCBLD2 Mouse Monoclonal Antibody
产品货号: AMM83487




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using DCBLD2 mouse mAb against PC-3 lysate



Lane 1: PC-3 lysate ; Lane 2: DCBLD2 immunoprecipitated from PC-3 lysate by AMM83487; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83487 can immunoprecipitate DCBLD2;

注意事项 (Note)

For research use only .