

产品名称: LXN Mouse Monoclonal Antibody
产品货号: AMM83498



产品概述 (Summary)

产品名称 (Production Name)	LXN Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LXN
别名 (Alternative Names)	LXN
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BS40.Purified recombinant fragment of human LXN

产品应用 (Application)

稀释比 (Dilution Ratio)	IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: LXN Mouse Monoclonal Antibody
产品货号: AMM83498



研究领域 (Research Area)

图片 (Image Data)



Lane 1: A549 lysate; Lane 2: LXN immunoprecipitated from A549 lysate by AMM83498; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83498 can immunoprecipitate LXN;

注意事项 (Note)

For research use only .