

产品名称: GSTT2B Mouse Monoclonal Antibody
产品货号: AMM83520



产品概述 (Summary)

产品名称 (Production Name)	GSTT2B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTT2B
别名 (Alternative Names)	GSTT2B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P0CG30.Purified recombinant fragment of human GSTT2B

产品应用 (Application)

稀释比 (Dilution Ratio)	IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: GSTT2B Mouse Monoclonal Antibody
产品货号: AMM83520



研究领域 (Research Area)

图片 (Image Data)



Lane 1: Rat liver tissue lysate; Lane 2: GSTT2B immunoprecipitated from rat liver tissue lysate by AMM83520; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83520 can immunoprecipitate GSTT2B;

注意事项 (Note)

For research use only .