

产品名称: TRIM49 Mouse Monoclonal Antibody
产品货号: AMM83531



产品概述 (Summary)

产品名称 (Production Name)	TRIM49 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIM49
别名 (Alternative Names)	TRIM49
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P0CI25.Purified recombinant fragment of human TRIM49

产品应用 (Application)

稀释比 (Dilution Ratio)	IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: TRIM49 Mouse Monoclonal Antibody
产品货号: AMM83531




研究领域 (Research Area)

图片 (Image Data)



Lane 1: HEK-293 lysate (transfected with DYKDDDDK tagged TRIM49 (Cat. P67967QATN) at N-terminus); Lane 2: TRIM49 immunoprecipitated from the same lysate as Lane 1; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83531 can immunoprecipitate TRIM49;



Lane 1: HEK-293 lysate (transfected with EGFP tagged TRIM49 (Cat. P67967QBTN) at C-terminus); Lane 2: TRIM49 immunoprecipitated from the same lysate as Lane 1; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83531 can immunoprecipitate TRIM49;

注意事项 (Note)

For research use only .