

产品名称: MYH9 Mouse Monoclonal Antibody
产品货号: AMM83577



产品概述 (Summary)

产品名称 (Production Name)	MYH9 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH9
别名 (Alternative Names)	MYH9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P35579.Purified recombinant fragment of human MYH9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:250,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	227kDa


研究背景 (Background)

产品名称: MYH9 Mouse Monoclonal Antibody
产品货号: AMM83577




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 μ g of HeLa lysate; Lane 2: 15 μ g of HT-29 lysate; Lane 3: 15 μ g of A-431 lysate; Result: AMM83577 can detect MYH9 by Western blotting;



Lane 1: HeLa lysate; Lane 2: MYH9 immunoprecipitated from HeLa lysate by AMM83577; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83577 can immunoprecipitate MYH9;



IHC-P analysis of kidney tissue by MYH9 antibody (AMM83577)

产品名称: MYH9 Mouse Monoclonal Antibody
产品货号: AMM83577



注意事项 (Note)

For research use only .