

产品名称: SMS Mouse Monoclonal Antibody
产品货号: AMM83579



产品概述 (Summary)

产品名称 (Production Name)	SMS Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMS
别名 (Alternative Names)	SMS
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P52788.Purified recombinant fragment of human SMS

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: SMS Mouse Monoclonal Antibody
产品货号: AMM83579




研究领域 (Research Area)

图片 (Image Data)



Lane 1: A-375 lysate; Lane 2: SMS immunoprecipitated from A-375 lysate by AMM83579; Lane 3: The same Lane 2 but use IgG isotype control antibody; Result: AMM83579 can immunoprecipitate SMS;



IHC-P analysis of testis tissue by SMS antibody (AMM83579)

注意事项 (Note)

For research use only .