

产品名称: Human IgE Mouse Monoclonal Antibody
产品货号: AMM83587



产品概述 (Summary)

产品名称 (Production Name)	Human IgE Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Human IgE
别名 (Alternative Names)	Human IgE
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01854.Recombinant human IgE

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	47kDa


研究背景 (Background)

产品名称: Human IgE Mouse Monoclonal Antibody
产品货号: AMM83587



研究领域 (Research Area)

图片 (Image Data)



100 ng of human IgE was run on 6-18% SDS-PAGE under reducing conditions and blotted onto nitrocellulose membrane. AMM83587 at 1 μ g/mL was used as the primary antibody and peroxidase conjugated goat anti-mouse IgG was used as the secondary antibody. Human IgE band was visualized using ECL Substrate.

注意事项 (Note)

For research use only .