

产品名称: ASRGL1 Mouse Monoclonal Antibody
产品货号: AMM83608



产品概述 (Summary)

产品名称 (Production Name)	ASRGL1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASRGL1
别名 (Alternative Names)	ASRGL1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q7L266.Purified recombinant fragment of human ASRGL1

产品应用 (Application)

稀释比 (Dilution Ratio)	IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: ASRGL1 Mouse Monoclonal Antibody
产品货号: AMM83608



研究领域 (Research Area)

图片 (Image Data)



Lane 1: Hep G2 lysate; Lane 2: ASRGL1 immunoprecipitated from Hep G2 lysate by AMM83608; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83608 can immunoprecipitate ASRGL1;

注意事项 (Note)

For research use only .