

产品名称: SLC4A1 Mouse Monoclonal Antibody
产品货号: AMM83609



产品概述 (Summary)

产品名称 (Production Name)	SLC4A1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC4A1
别名 (Alternative Names)	SLC4A1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P02730.Purified recombinant fragment of human SLC4A1

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: SLC4A1 Mouse Monoclonal Antibody
产品货号: AMM83609




研究领域 (Research Area)

图片 (Image Data)



Lane 1: Kidney tissue lysate ; Lane 2: SLC4A1 immunoprecipitated from kidney tissue lysate by AMM83609 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83609 can immunoprecipitate SLC4A1;



IHC-P analysis of kidney tissue by SLC4A1 antibody (AMM83609)

注意事项 (Note)

For research use only .