

产品名称: KHK Mouse Monoclonal Antibody
产品货号: AMM83610



产品概述 (Summary)

产品名称 (Production Name)	KHK Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KHK
别名 (Alternative Names)	KHK
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50053.Purified recombinant fragment of human KHK

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: KHK Mouse Monoclonal Antibody
产品货号: AMM83610




研究领域 (Research Area)

图片 (Image Data)



Lane 1: Rat kidney tissue lysate; Lane 2: KHK immunoprecipitated from rat kidney tissue lysate by AMM83610; Lane 3: The same Lane 2 but use IgG isotype control antibody; Result: AMM83610 can immunoprecipitate KHK;



IHC-P analysis of liver tissue by KHK antibody (AMM83610)

注意事项 (Note)

For research use only .