

产品名称: DMXL2 Mouse Monoclonal Antibody
产品货号: AMM83617



产品概述 (Summary)

产品名称 (Production Name)	DMXL2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DMXL2
别名 (Alternative Names)	DMXL2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8TDJ6.Purified recombinant fragment of human DMXL2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	339.6kDa


研究背景 (Background)

产品名称: DMXL2 Mouse Monoclonal Antibody
产品货号: AMM83617




研究领域 (Research Area)

图片 (Image Data)



Lane 1: 15 μ g of Jurkat lysate; Lane 2: 15 μ g of HeLa lysate; Lane 3: 15 μ g of HEK-293 lysate; Result: AMM83617 can detect DMXL2 by Western blotting;



Lane 1: HEK-293 lysate; Lane 2: DMXL2 immunoprecipitated from HEK-293 lysate by AMM83617; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83617 can immunoprecipitate DMXL2;

注意事项 (Note)

For research use only .