

产品名称: TSNAX Mouse Monoclonal Antibody  
产品货号: AMM83620



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | TSNAX Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody       |
| 宿主 (Host)              | Mouse                           |
| 应用 (Application)       | WB,IP                           |
| 种属反应性 (Reactivity)     | Human                           |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                     |
|------------------------|-----------------------------------------------------|
| 基因名 (Gene Name)        | TSNAX                                               |
| 别名 (Alternative Names) | TSNAX                                               |
| 基因 ID (Gene ID)        |                                                     |
| 蛋白 ID (SwissProt ID)   | Q99598.Purified recombinant fragment of human TSNAX |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IP 1:100-1:200             |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 33 kDa; Observed MW: 33 kDa |


## 研究背景 (Background)

产品名称: TSNAX Mouse Monoclonal Antibody  
产品货号: AMM83620




## 研究领域 (Research Area)

## 图片 (Image Data)



Lane 1: 15  $\mu$ g of MCF7 lysate; Lane 2: 15  $\mu$ g of HeLa lysate; Lane 3: 15  $\mu$ g of HEK-293 lysate; Lane 4: 15  $\mu$ g of SK-OV-3 lysate;

Result: AMM83620 can detect TSNAX by Western blotting.



Lane 1: TSNAX immunoprecipitated from MCF7 lysate by P154562; Lane 2: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83620 can immunoprecipitate TSNAX.

## 注意事项 (Note)

For research use only .